

Somente estrutura

| Altura A (mm) | Largura L (mm) | Profundidade P (mm) | | | |
|---------------|----------------|---------------------|-------------------|-------------------|-------------------|
| | | 400 | 600 | 800 | 1000 |
| 1500 | 610 | GP 15064.E | GP 15066.E | GP 15068.E | |
| 1700 | 610 | GP 17064.E | GP 17066.E | GP 17068.E | |
| 1900 | 610 | GP 19064.E | GP 19066.E | GP 19068.E | GP 19060.E |
| 2100 | 610 | GP 21064.E | GP 21066.E | GP 21068.E | GP 21060.E |
| 2200 | 610 | GP 22064.E | GP 22066.E | GP 22068.E | GP 22060.E |
| | | | | | |
| 1500 | 810 | GP 15084.E | GP 15086.E | GP 15088.E | |
| 1700 | 810 | GP 17084.E | GP 17086.E | GP 17088.E | |
| 1900 | 810 | GP 19084.E | GP 19086.E | GP 19088.E | GP 19080.E |
| 2100 | 810 | GP 21084.E | GP 21086.E | GP 21088.E | GP 21080.E |
| 2200 | 810 | GP 22084.E | GP 22086.E | GP 22088.E | GP 22080.E |
| | | | | | |
| 1700 | 1010 | GP 17104.E | GP 17106.E | GP 17108.E | |
| 1900 | 1010 | GP 19104.E | GP 19106.E | GP 19108.E | GP 19100.E |
| 2100 | 1010 | GP 21104.E | GP 21106.E | GP 21108.E | GP 21100.E |
| 2200 | 1010 | GP 22104.E | GP 22106.E | GP 22108.E | GP 22100.E |
| | | | | | |
| 1700 | 1210 | GP 17124.E | GP 17126.E | GP 17128.E | |
| 1900 | 1210 | GP 19124.E | GP 19126.E | GP 19128.E | GP 19120.E |
| 2100 | 1210 | GP 21124.E | GP 21126.E | GP 21128.E | GP 21120.E |
| 2200 | 1210 | GP 22124.E | GP 22126.E | GP 22128.E | GP 22120.E |

Coluna Auxiliar

| Altura A (mm) | Largura L (mm) | Profundidade P (mm) | | | |
|---------------|----------------|---------------------|-------------------|-------------------|-------------------|
| | | 400 | 600 | 800 | 1000 |
| 1900 | 410 | GP 19044.0 | GP 19046.0 | GP 19048.0 | GP 19040.0 |
| 2100 | 410 | GP 21044.0 | GP 21046.0 | GP 21048.0 | GP 21040.0 |
| 2200 | 410 | GP 22044.0 | GP 22046.0 | GP 22048.0 | GP 22040.0 |

A= Altura do gabinete sem a base.